PRODUCTS ARE SOLD FOR INSTRUCTIONAL USE ONLY AND ARE NOT FOR GENERAL CONSUMER, PERSONAL, DIRECT OR INDIRECT HUMAN OR ANIMAL USE.

IDENTITY (As Used on Label and List) Demo Dose[®] Simulated Blood Powder Gallon and Quart

_			
	Section	I	

Manufacturer's Name Pocket Nurse [®] Emergency Telephone Number Chemtrec (USA) 800-424-9300 Canutec (Canada) 613-996-6666 Outside USA & Canada 703-527-3887

Address (Number, Street, City, State, and ZIP Code) 200 1st Street Ambridge Pa 15003 Telephone Number for Information 800-225-1600 Date Prepared February 01, 2006 Signature of Preparer (Optional)

Hazardous Components (Specific Chemical Identity; Common Name(s))	CAS #	OSHA PEL	ACGIH TLV	Other Limits	% Optional
Inorganic Salts – Non Hazardous (proprietary information)					5 %
FD&C Red #40 Gran Color Additive – Non Hazardous	25956-17-6				95 %

Section II - Hazardous Ingredients/Identity Information

Section III - Physical/Chemical Characteristics

Boiling Point	N/A	рН	N/A
Vapor Pressure (mm Hg.)	N/A	Melting Point	N/A
Vapor Density (AIR = 1)	N/A	Evaporation Rate n-butyl acetate = 1	N/A

Solubility in Water: Completely Soluble

Appearance and Odor: Dry Colorant; odorless

Section IV - Fire and Explosion Hazard

Flash Point (Method Used)	Flammable Limits	LEL	UEL
N/A	N/A	N/A	N/A

Extinguishing Media

Water, chemical foam or dry chemical extinguisher

Special Fire Fighting Procedures None

Unusual Fire and Explosive Hazards

None

Section V - Reactivity Data

Stability	Unstable		Conditions to Avoid: None
	Stable	Х	

Incompatibility (Materials to Avoid) None

Hazardous Decomposition or Byproducts None

Hazardous Polymerization	May Occur		Conditions to Avoid
	Will Not Occur	Х	N/A

Section VI - Health Hazard Data

Ingestion: None Eye Contact: None Skin Contact: None Inhalation: None Other: None

Signs and Symptoms of Exposure

Irritation of eyes may be possible with contact.

Emergency and First Aid Procedures

Eyes: Flush thoroughly with water, including upper and lower lids. If irritation develops or persists, get medical attention.

Section VII - Precautions for Safe Handling and Use

Steps to be taken in Case Material is Released or Spilled

If dry, vacuum or sweep color up for disposal. If wet, Absorb with paper towels, or other absorbent material. Place in suitable container. Slippery when wet.

Waste Disposal Method

Disposal may be subject to Federal, State/Provincial and local regulations. Check your local ordinances.

Precautions to be Taken in Handling and Storage

Store in a tightly sealed container. Store in cool, dry place.

Other Precautions

For educational use only. Not for clinical or diagnostic use. To be used only with direct supervision of a qualified professional specifically trained in laboratory procedures. This is not real Blood. Do not take or use internally. Do not inject. Wash hands thoroughly with soap and water after used. Not for drug, food or household use. Keep out of reach of children.

Section VIII - Control Measures

Respiratory Protection (Specify Type)

No respiratory protection required. This product can generate dust, so dust masks are suggested.

Ventilation	Local Exhaust none needed	Special None required	
	Mechanical None needed	Other None	
Protective Gloves	Nitrile	Eye Protection Eye protection is suggested for compliance with Good Manufacturing Practices	
		(GMP)	

Other Protective Clothing or Equipment

None Required. Eyewash station in proximity, within 15 sec. of work station is suggested.

Work/Hygienic Practices

Wear protective equipment. Wash hands thoroughly after use. May temporarily stain skin. Use good personal hygiene practices

NOTICE: Pocket Nurse® expressively disclaims all express or implied warranties of merchantability and fitness for a particular purpose, with respect to the product or information provided herein.

All information appearing herein is based upon data obtained from the manufacturer and/or recognized technical sources. While the information is believed to be accurate, Pocket Nurse® makes no representations as to its accuracy or sufficiency. Conditions of use are beyond Pocket Nurse® control and therefore users are responsible to verify this data under their own operating conditions to determine whether the product is suitable for their particular purposes and they assume all risks of their use, handling and disposal of the product or from the publication or use of, or reliance upon, information contained herein. This information relates only to the product designated herein, and does not relate to its use in combination with any other material or in any other process.